INNO-SEPT FRESH Disinfecting Solution

Microbial spectrum: bactericide (MRSA), yeasticide, virucide, tuberculocide.

Requirement for application: undiluted state.

Areas and methods of usage: the preparation can be used in various areas of public health care, food industry, catering and other services. It is designed for professional application.

Hygienic disinfection of hands: dose approximately 3 ml of the solution INNO-SEPT FRESH to the palm and rub it in thoroughly for 30 seconds (1 minute for tuberculocidal effect).

Disinfecting surfaces: spray the surface with the preparation or treat it with the wet wipe imbibed by the disinfectant, leaving no part of the surface without treatment. Minimum working time: 5 minutes. After the working time has elapsed, the treated surface can be rinsed with water of drinking-water quality. In case of stronger contaminations, we recommend the preliminary application of a

more intensive preparation for washing hands (MANUELA) and cleaning surfaces (INNOFLUID-ALK-M), respectively.





EU hazard symbol and letter designator:

Xi irritant,

N dangerous for the environment

Flammability degree: III

Safety instructions: do not apply the disinfecting solution INNO-SEPT FRESH for treating mucous membranes, running wounds or the skin surface near the eyes. Use the preparation exclusively for its function in the appropriate way. Comply with general labour safety and labour hygienic instructions.

R 10 - Flammable.

R 36 - Irritating to eyes.

R 51/53 - Toxic to aquatic organisms; may cause long-term adverse effects in the aquatic environment.

S 16 - Keep away from sources of ignition - no smoking.

S 25 - Avoid contact with eyes.

S 26 - In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S 61 - Avoid release to the environment - refer to special instructions / safety data sheet.

S 64 - If swallowed, rinse mouth with water (only if the person is conscious).

In case of accidental swallow, remove the victim to fresh air, give her/him plenty of water to drink and call a physician to the premises immediately. Do not induce vomiting!

Storage: store the preparation in its original unopened package at a dry cool but freeze-free place protected from sunshine, separated from foodstuffs, sources of ignition and heat.

S 2 - Keep out of the reach of children.

Environmental protection: do not discharge the disinfecting solution or its container to the environment. Handle the generated sewage according to local instructions, considering the quality of the contaminations.

Quality preservation time: 2 years after the production date in case of proper storage conditions.

Packaging units: 201 and 51 can, 11 flacon.

Produced and distributed by: INNOVENG 1 Kft. **Address:** H-1114 Budapest, Szabolcska M. u. 5.

Internet site: www.innoveng1.hu

Commercial office - address: H-2040 Budaörs, Liget utca 3/2; telephone: +36-23/80-17-80, +36-23/80-17-81; fax:

+36-23/421-933, +36-70/362-8756